

**IMCO Corporation Ltd AB**  
PO Box 21195  
SE -100 31 Stockholm  
Sweden  
[www.imcocorp.se](http://www.imcocorp.se)

## PRODUCT SPECIFICATION

<b>Product No.</b>	<b>AGRX-04</b>
<b>Product description:</b>	<b>Antibody (goat) to Human glutaredoxin 1, biotinylated</b>
<b>Lot No.</b>	<b>1</b>
<b>Specifications:</b>	<p>The gamma globulin fraction from goat, immunized with recombinant human glutaredoxin 1 (Grx1) was affinity purified on a column of human Grx1 and biotinylated.</p> <p>Each ampoule (192 µl) contains 50µg of antibody in PBS + 15 mM sodium azide, pH 7.5. The preparation (0.26 mg/ml) comprises inhibitory as well as non-inhibitory antibodies and is useful for ELISA applications. See protocol enclosed</p>
<b>Reference:</b>	Lundberg, M., Potamitou-Fernandes, A., Kumar, S. and Holmgren, A.: Cellular and plasma levels of human glutaredoxin 1 and 2 detected by sensitive ELISA systems. <i>Biochem. Biophys. Res. Commun.</i> , <u>319</u> , 801-809, 2004
<b>Storage:</b>	Keep at +4°C
<b>Examples of use:</b>	ELISA see protocol Good for Western blotting . Good for Immunohistochemistry.

**NOTE!** During shipment, the contents of the antibody ampoule may have been redistributed on the inner walls. Before opening the screw-cap of the antibody tube, we recommend to centrifuge the tube shortly, in order to recover liquid from tube cap and walls.

## Grx 1 ELISA Protocol

### Buffers:

Sodium Carbonate buffer pH 9.7

0.1% Tween in PBS (PBS-T)

PBS, 0.5% BSA, 0.02% NaN<sub>3</sub> (incubation buffer)

10% Diethanolamine, 0.5mM MgCl<sub>2</sub>, 0.02% NaN<sub>3</sub> pH 9.8 (substrate buffer)

### Day 1

1. Dilute anti-Grx1 (goat) to 5 µg/ml in Sodium Carbonate buffer pH 9.7.
2. Coat the wells by adding 50 µl 5 µg/ml anti-Grx1 into the wells.
3. Incubate the plate o/n at 4 °C with slight shake.
4. Prepare sample and Grx1 standard

#### *Grx1 Standard:*

Prepare 200 µl 50 ng/ml Grx1 in incubation buffer.

#### *Sample:*

Mix 100 µl sample with 100 µl incubation buffer.

5. Incubate the plate o/n at 4 °C with slight shake.

### Day 2

6. Discard the content of the plate.
7. Add 200 µl in the wells and incubate 3 hours at 4 °C with slight shake.
8. Discard the incubation buffer and wash the plate 4x with PBS-T.
9. Add 50 µl incubation buffer in those wells where standards and samples are to be diluted.
10. For standards; add 100 µl in wells A12 and B12, remove 50 µl and add it to the next well and so on.
11. For the sample; add 100 µl in the desired wells, remove 50 µl and add it to the next well and so on.
12. Incubate the plate o/n at 4 °C with slight shake.

### Day 3

13. Discard the content of the plate and wash the plate 4x with PBS-T.
14. Dilute biotinylated anti-Grx1 (goat) to 1 µg/ml in incubation buffer and add it to the wells.
15. Incubate 2 hours at room temperature with slight shake.
16. Discard the content of the plate and wash the plate 4x with PBS-T.
17. Dilute streptavidin-ALP 1:1000 in incubation buffer, add 50 µl to each well incubate for 1 hour at room temperature with slight shake.
18. While incubating; add 1 tablet 4-nitrophenyl phosphate in 5 ml substrate buffer.
19. Discard the content of the plate and wash the plate 6x with PBS-T.
20. Add 50 µl 4-nitrophenyl phosphate to the wells.
21. Incubate at least 15 min (sometimes the reaction takes longer)
22. Read A<sub>405</sub>.