

IMCO Corporation Ltd AB

PO Box 21195
SE -100 31 Stockholm
Sweden
www.imcocorp.se

PRODUCT SPECIFICATION

Product No.	FkGRX-01FS
Product description:	Fluorescent substrate for glutaredoxin assay with FkGRX-01
Lot No.	02
Specifications:	Fluorescent substrate Lyophilized fluorescent substrate <i>Preparation procedure:</i> Dissolve the content in 1.0 ml distilled water. Store at -20°C.
Reference:	See methods description in www.imcocorp.se
Storage:	Keep at -20°C
Examples of use:	A substrate for measurement of glutaredoxin using the FkGRX-01 kit from IMCO. See www.imcocorp.se

NOTE! Make sure to reconstitute all of the lyophilized material in the ampoule. After adding water, close the screw-cap again, shake vigorously, then centrifuge shortly, in order to recover liquid from tube cap and walls.