

IMCO Corporation Ltd AB

PO Box 21195
SE -100 31 Stockholm
Sweden
www.imcocorp.se

PRODUCT SPECIFICATION

Product No.	FkTRX-02-V2-FS
Product description:	Fluorescent substrate for Trx assay with FkTRX-02-V2
Lot No.	01
Specifications:	Fluorescent substrate Lyophilized fluorescent substrate <i>Preparation procedure:</i> Dissolve the content in 1.0 ml distilled water. Store at -20°C
Reference:	See methods description for FkTRX-02-V2 at www.imcocorp.se
Storage:	Keep at -20°C
Examples of use:	A substrate for measurement of thioredoxin activity using FkTRX-02-V2 kit from IMCO. See www.imcocorp.se

NOTE! Make sure to reconstitute all of the lyophilized material in the ampoule. After adding water, close the screw-cap again, shake vigorously, then centrifuge shortly, in order to recover liquid from tube cap and walls.